

**Bioinformatics and Bioengineering Summer Institute (2003)**  
**Introduction to the Institute – Dilution Problems**

1. You want to dilute a culture of *E. coli* at a concentration of  $2 \cdot 10^9$  cells/ml to a point where you can spread 100  $\mu$ l on a plate and get a countable number of colonies (define "countable" as 50 to 300).
  - 1a. What volume of the original culture would you need to pipette into 100  $\mu$ l fresh growth medium in order to get a countable number of colonies?
  - 1b. Considering that the smallest volume you can pipet accurately is about 2  $\mu$ l, can you perform the dilution calculated in 4a?
  - 1c. Fill in the blanks in the dilution scheme below, designed to achieve the desired objective:

Culture	Operation	Concentration	Cells in 100 $\mu$ l (= cells on plate)
Culture 0	(original culture)	$2 \cdot 10^9$ cells/ml	$2 \cdot 10^8$ cells
Culture A	1 ml medium + 10 $\mu$ l of Culture 0	$2 \cdot 10^7$ cells/ml	$2 \cdot 10^6$ cells
Culture B	1 ml medium + 10 $\mu$ l of Culture A		
Culture C	1 ml medium + 5 $\mu$ l of Culture B		

2. You want to dilute a culture of *E. coli* at a concentration of  $1 \cdot 10^9$  cells/ml to a point where you can spread 200  $\mu$ l on a plate and get a countable number of colonies (define "countable" as 50 to 300).
  - 2a. What volume of the original culture would you need to pipette into 100  $\mu$ l fresh growth medium in order to get a countable number of colonies?
  - 2b. Considering that the smallest volume you can pipet accurately is about 2  $\mu$ l, can you perform the dilution calculated in 4a?
  - 2c. Fill in the blanks in the dilution scheme below, designed to achieve the desired objective:

Culture	Operation	Concentration	Cells in 100 $\mu$ l (= cells on plate)
Culture 0	(original culture)	$1 \cdot 10^9$ cells/ml	$2 \cdot 10^8$ cells
Culture A	1 ml medium + 10 $\mu$ l of Culture 0	$1 \cdot 10^7$ cells/ml	$2 \cdot 10^6$ cells
Culture B	1 ml medium + 10 $\mu$ l of Culture A		
Culture C	1 ml medium + 5 $\mu$ l of Culture B		

3. You wish to grow *E. coli* in 5 ml of liquid medium containing 15  $\mu\text{g/ml}$  of the antibiotic tetracycline. You have 5 ml of liquid medium without antibiotic and you have a stock solution of tetracycline at a concentration of 15 mg/ml.
  - 3a. By what factor must the 15 mg/ml stock solution of tetracycline be diluted to reach the final concentration of 15  $\mu\text{g/ml}$ ? Call that factor the "dilution factor".
  - 3b. What volume of the stock solution must you add to 5 ml in order for the tetracycline to be diluted by the dilution factor?
  
4. You wish to grow *E. coli* in 10 ml of liquid medium containing 50  $\mu\text{g/ml}$  of the antibiotic kanamycin. You have 10 ml of liquid medium without antibiotic and you have a stock solution of kanamycin at a concentration of 25 mg/ml.
  - 4a. By what factor must the 25 mg/ml stock solution of kanamycin be diluted to reach the final concentration of 50  $\mu\text{g/ml}$ ? Call that factor the "dilution factor".
  - 4b. What volume of the stock solution must you add to 10 ml in order for the kanamycin to be diluted by the dilution factor?
  
5. **Person #1:** Using a micropipettor, pipet 150  $\mu\text{l}$  of water 3 times to a weighing dish.  
**Person #2:** Using the same micropipettor, pipet 75  $\mu\text{l}$  of water 6 times from the weighing dish. (you should exactly exhaust the water).
  
6. **Person #1 and #2:** Using a glass pipet, suck up 6.6 ml water and distribute 2.2 ml to each of three weighing dishes.